Instruction For Use LP2-A060 Version G

Accuracy

Real Time

Radiation
Free

Safety

Educational
Tool

Table of content

EN - Instructions For Use, Assembly Instructions and Surgical Technique of the DSG® Integration Module for Screw

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14

In no case can the DSG® Integration Module for Screw replace the surgeon's experience or knowledge of anatomic structures. It is intended to assist the surgeon in the decision-making process during surgery in the operating room.

This User Manual is a combination of the *Instructions for Use* document, the *Surgical Technique* document and the *Assembly Instruction* document. All of the essential information from each of those source documents has been faithfully presented in this User Manual.

1. INTRODUCTION

The DSG® Integration Module for Screw, combined with devices from the established partner(s) of SpineGuard, is used in drilling pilot holes into the spinal vertebrae. This system serves to alert the surgeon to a possible vertebral cortex fracture while inserting a pedicle screw.

The DSG® integration module for screw assembly is the combination of the DSG® pin connected to the Electronic Handle attached to the instrument from the established partner(s) through the DSG® ratcheting handle.

It is strongly recommended that the working principle of the DSG® Technology be clearly understood before using the DSG® Integration Module for Screw.

2. TECHNICAL DESCRIPTION

The DSG® Integration Module for Screw is part of the DSG® product range. The DSG® Integration Module for Screw consists of:

- DSG® single-use instruments provided sterile:
 - DSG® Pin Ref. D1PUXXXX, which embeds the proprietary DSG® bipolar sensor
 - Electronic Handle Ref. P2HEXXXX, which houses the electronics that read and translate the electrical signals detected by the sensor.
Note: An optional feature is available to transmit the conductivity data to an external display to allow signal visualization and recording. The handle equipped with this feature is labelled DSG Connect and is equipped with a communication module which transmits data through radio frequency (RF) communication (2.40 – 2.48 GHz, GFSK modulation and maximum RF output power of 9.9 dBm)
- DSG® reusable instruments that must be cleaned and sterilized before use:
 - DSG® Ratcheting Handle - Ref. D1RHXXXX
 - DSG® Technology Instrument Tray: Base Ref. D1TT1XXX and Lid Ref. D1TT2XXX.
- Optional DSG Connect visualization software: DSG Connect App – Ref D1SW0001 to be used on a dedicated tablet and with the Electronic handle equipped with the DSG Connect feature

The DSG® Integration Module for Screw is meant to be used with a cannulated pedicle screw and a screwdriver from an established partner.

The tip of the instrument is the applied part.

Note: the DSG® Integration Module for Screw is a battery-operated device with a maximum output current of 5.5 milliamps ("mA"). Published literature suggests that threshold level monitoring requires at least 8-10 mA. The DSG® Integration Module for Screw, which does not perform threshold recordings, is not intended to be used as a substitute for a threshold level monitoring device.

Essential Performances:

- Mechanics essential performance: the cutting ability during the pedicle drilling.
- Electronics essential performance: sound and visual signals lead the surgeon during the pedicle drilling. Those informing signals convey the conductivity information which characterizes the tissues in contact with the sensor tip.

3. BEFORE USE

This is the inspection, cleaning and sterilization information for the DSG® Integration Module for Screw.

3.1 Information related to single-use devices

The electrical and sensory components of the DSG® Integration Module for Screw are single use devices. The retreatment and the re-use of the single-use DSG® Integration Module for Screw components are forbidden.

Inspect the single-use instruments for integrity of packaging. ONLY USE INSTRUMENTS WITH INTACT PACKAGING.

DO NOT CLEAN, RE-STERILIZE, OR RE-USE THE SINGLE USE COMPONENTS OF THIS DEVICE.

Several potential dangers related to the re-use of the DSG® Integration Module for Screw can increase the risk to the patient. Some technical characteristics of the electrical and sensory components of the DSG® Integration Module for Screw are not compatible with the re-use of the instrument. Possible issues with re-sterilization and re-use are: non-functioning device, failure of some or all of the electronic functions, unknown changes to the electronic function of the device, non-sterile device, transfer of disease, cleaning agents or disinfectants.

Additionally, the re-use of a single use device will result in the loss of traceability of the medical device and the loss of the technical documentation, such as this instruction for use.

3.2 Information related to reusable devices

Preliminary Check

Before each use, all DSG® reusable instruments must be inspected, cleaned and sterilized.

All assembled instruments must be disassembled before inspection, cleaning and sterilization. To disassemble:

1. If still present, remove the Electronic Handle and the DSG® Pin from the DSG® Ratcheting Handle. Dispose of the Electronic Handle and the DSG® Pin (single use only).
2. If still assembled, separate the partner screwdriver from the DSG® Ratcheting Handle.

Prior to each use, perform a visual inspection of all the disassembled reusable instruments.

Ensure that the V-Shape in the DSG® Ratcheting Handle is not damaged (twisted or broken). Refer to Fig. 1


Fig.1

If any damage or excessive wear is noted on the DSG® reusable instruments, such as twisting, bending, breakage, pinching at tips, lack of integrity of the etch marks, discoloration or cracks, they should not be used and should be discarded following the hospital procedures for post-treatment of contaminated devices.

The ratcheting handle should be lubricated after each use with 1 drop of standard biocompatible autoclavable oil at the marked points (see fig.2). Lubricate the locking ring and the selector of the ratcheting mechanism and distribute the oil by moving the components back and forth several times. Wipe off the excess oil with a cloth.


Fig.2

Cleaning for the DSG® reusable instruments

Disassembly of the instruments for cleaning should be performed in the designated cleaning area immediately following use at the point of care. The DSG® reusable instruments are delivered clean but not sterile; therefore, they must be sterilized before use. The surfaces of all DSG® instruments must be cleaned and decontaminated before sterilization. Clean all DSG® reusable instruments after each surgery at the point of use with a disposable wipe to remove excessive soil.

For the initial and any subsequent cleaning, the recommended general protocol for the pre-disinfection, decontamination, and cleaning of all these instruments is as follows:

- In a pre-disinfection bath, immerse and soak the instruments in an enzymatic detergent solution (pH<11, equivalent to «Steris® Prolystica 2X Enzymatic cleaner») prepared in accordance to the manufacturer's instructions (0.2 to 0.8%) for a minimum of 15 minutes at room temperature (15-25°C).
- After soaking, use a soft-bristled nylon brush to gently scrub the instruments for at least 1 minute, until all visible soil has been removed. Actuate any moveable mechanisms while cleaning. Pay particular attention to lumen and all difficult-to-reach areas. To thoroughly scrub the lumen, push a long, narrow, soft nylon brush in and out of the lumen using a twisting motion.
- Use a syringe, pipette or water pistol to thoroughly flush lumen and all difficult-to-reach areas with an enzymatic detergent solution (pH<11, prepared in accordance to the manufacturer's instructions).
- Ultrasonically clean the instruments (completely immersed) for 15 minutes in enzymatic detergent solution (pH<11, prepared in accordance to the manufacturer's instructions) at ambient temperature (15-25°C).
- Rinse the instruments in distilled water for at least 1 minute at room temperature (15- 25°C). Actuate any moveable parts while rinsing. Thoroughly flush lumens and all difficult-to-reach areas.
- Visually inspect the device and repeat the manual pre-cleaning if there is any visible residual debris.
- Load the instruments into a validated washer/disinfector so that lumen can drain.
- Run an automatic cleaning cycle using the minimum parameters below:

Cycle	Minimum Time	Minimum Temperature	Type of Detergent / Water
Pre-cleaning	2 minutes	Cold < 45°C	Tap Water
Cleaning	10 minutes	Heated 55°C	Enzymatic Detergent (pH< 11, equivalent to «Steris® Prolystica 2X Enzymatic cleaner» 0.2 to 0.8%)
Rinse 1	2 minutes	Cold < 45°C	Tap Water
Rinse 2	2 minutes	Cold < 45°C	Tap Water
Thermal Rinse	5 minutes	Heated 90°C	High Purity Water*
Dry	20 minutes	Heated > 60°C	Not applicable

*Water extensively treated (usually by a multistep treatment process that may include a carbon bed, softening, DI, and RO or distillation) to ensure that the microorganisms and the inorganic and organic material are removed from the water (AAMI TIR 34)

- When unloading, visually inspect all instruments. Repeat the manual and automated washing if there is any residual debris.

Note: Certain washing solutions such as those containing bleach or formalin may damage the instruments and should not be used.

Note: Since no reprocessing methods have been validated for removal or inactivation of transmissible spongiform encephalopathy (TSE) agents (i.e., prions) from medical devices, this device should not be used for patients with known or suspected TSE agent disease.

The cleaning instructions provided by SpineGuard® were validated to demonstrate the adequacy of the cleaning instructions and methods. Users are separately responsible for implementing reprocessing following the instructions, training employees, and monitoring their ability to reprocess devices.

Sterilization for the DSG® reusable instruments

Before use, the disassembled and cleaned reusable components of the DSG® Integration Module for Screw must be placed in the DSG® Technology Instrument Tray for steam sterilization with the maximum following quantity:

Maximum quantity of instruments per DSG® Technology Instrument Tray	
DSG® Ratcheting Handles	2

The Tray should be wrapped with a double layer of sterilization wrap. Use a validated, properly maintained and calibrated steam sterilizer. Effective steam sterilization can be achieved using the following cycle:

GEOGRAPHIC AREA	METHOD	CYCLE TYPE	TEMPERATURE	EXPOSURE TIME	DRY TIME
Outside the U.S.A	Steam	Pre-Vacuum	274°F (134 °C)	18 minutes	>15 minutes

THESE REPROCESSING INSTRUCTIONS HAVE BEEN VALIDATED AS BEING CAPABLE OF PREPARING THE DSG® REUSABLE INSTRUMENTS FOR REUSE. THEREFORE, SPINEGUARD® ASSUMES NO RESPONSIBILITY FOR INAPPROPRIATE REPROCESSING BY THE USER.

4. USE

4.1 Selection of the appropriate DSG® Integration Module for Screw components

The DSG® Integration Module for Screw reusable instruments must be used only with the DSG® Pin, the Electronic Handle and the instruments provided by the partner(s). Select the adequate DSG® Pin, and DSG® Ratcheting Handle according to the partner. The compatibility information is provided on the DSG® Pin packaging.

4.2 Assembly Instructions and functional check

Set-up instructions:

1. Attach the selected partner screwdriver to the Ratcheting Handle (Fig.3).

Fig.3

2. Check that the instrument is fully engaged with the Ratcheting Handle – the alignment mark on the partner screwdriver should be flush with the distal end of the Ratcheting Handle. This ensures the complete engagement of the instrument with the Ratcheting Handle.


Fig.4

3. Insert the Pin through the Ratcheting Handle – partner screwdriver assembly (from step 1)

Fig.5

4. Check that the knob of the Pin is fully threaded into the Ratcheting Handle.
5. Connect the selected screw onto the screwdriver/ DSG® Integration Module for Screw assembly. The distal tip of the Pin should be visible (3mm +/- 1mm) at the distal end of the screw.


Fig.6

Ensure that the ratchet mechanism is fully functional after the Pin has been completely threaded. If necessary, adjust engagement of Pin with Ratcheting Handle.

6. Activate the DSG® Integration Module for Screw by pulling the contact tab out of the Electronic Handle. Activation is confirmed by a beeping sound and green LEDs lighting up. Connect the Electronic Handle to the Pin – DSG® screw module assembly. Ensure that the Electronic Handle is securely connected to the ratcheting handle (Fig.7-8). Once ON, DSG® Integration Module for Screw cannot be turned off. The DSG® Integration Module for Screw battery provides sufficient autonomy to perform a complete spinal surgery.

Fig.7

Fig.8

Note: The devices equipped with the DSG Connect feature can be paired to an external display. The pairing instructions and specific instructions for use are available within the DSG Connect App.

7. Check the functionality of the DSG® Integration Module for Screw assembly by placing the distal tip in saline solution (0.9% sodium chloride). A high pitch and high cadence sound should be heard. The green LEDs should flash at the same frequency. This confirms that the DSG® Integration Module for Screw is fully functional and ready for use. If a metal bowl is used, do not touch the bowl with the tip of the instrument. If no sound is heard, it may be necessary to clean the tip of the instrument using a mild scrubber (e.g., electrocautery scrubber).
8. The DSG® Integration Module for Screw is now ready for use

The DSG® Integration Module for Screw has several safety features which allow easy checking of proper operation:

Green LEDs	ON	Flashing	OFF	OFF	OFF
Yellow LEDs	OFF	OFF	OFF	OFF	ON
Sound	None	Rated sound	ON	None	None
Case	①	②	③	④	⑤

- ① The instrument is ON. Either the tip of the instrument does not make contact with tissue or is in contact with bone that has very low conductivity.
- ② The tip of the instrument is in contact with tissue.
- ③ The battery is too low: do not use the instrument.
- ④ Either the instrument is not ON, or the battery does not work. Verify that the instrument is ON. If the instrument is ON but the LEDs remain OFF, do not use the instrument.

- ⑤ The instrument tip is not functioning properly: do not use the instrument.

If the LEDs are not visible through the clear wall of the Electronic Handle, do not use the DSG® Integration Module for Screw.

4.3 Surgical Technique – Open and Mini-Open Techniques

Screw length and diameter selection and connection onto screwdriver: refer to recommendations provided by the partner

The surgeon should prepare the intended entry points for the pedicle screws using preferred open or mini open techniques to expose the posterior surface of the spine. The correct location for the entry point and angle of entry may be checked by using fluoroscopy.

In the case of a mini-open procedure, use the retractor provided with the pedicle screw instrumentation set until the entry point can be directly visualized. Prior to inserting the screw in the pedicle with the DSG® Integration Module for Screw, prepare the entry point of the pedicle and perforate the dorsal cortex using standard techniques.

Before and during initial insertion of the DSG® Integration Module for Screw into the pedicle, ensure that the ratcheting mechanism is in the neutral position. This will enhance the ease of penetration into the pedicle (Fig.11).

Fig.11

Direct the tip of the DSG® Integration Module for Screw to the soft tissue (wound) to initiate a rated high-pitched “beep”.

Place the tip of the DSG® Integration Module for Screw at the entry point of the pedicle.

As the DSG® Integration Module for Screw is inserted past the entry point into the pedicle, a high pitch and high cadence feedback may be heard. This is due to the tip being advanced through the blood collected at the entry point (Fig.12-13).

Fig.12


Fig.13

Pay close attention to the audio feedback provided as the DSG® Integration Module for Screw is being advanced while applying constant pressure. Once a few screw threads have come in contact with bone, engage the ratchet mechanism by turning it clockwise (Fig.14). This will allow to advance the instrument down the pedicle in small increments.


Fig.14

Advance the DSG® Integration Module for Screw along the targeted screw trajectory. As the device is being advanced into the pedicle through cancellous bone, a medium pitch and medium cadence audio signal must be heard all along the screw insertion.

If there is a decrease in pitch and cadence or a significant increase before the screw is fully inserted it is indicating that the tip of the screw is approaching the cortical wall with a risk of imminent cortical breach. In that case, the screw should not be further inserted and redirected. If the screw trajectory is satisfying but the screw is estimated too long, it should be replaced by a shorter screw.

Muscular contractions may occur if the cortex has been fractured. In this case, the surgeon must stop advancement of the screw, and check and/or alter the initial trajectory.

After each screw insertion, clean the tip of the DSG® Pin with a soft cloth impregnated with saline solution.

In Case Redirection Is Necessary

A change in pitch and cadence of the audio feedback indicates a change in the type of tissue around the tip of the DSG® Integration Module for Screw. If the audio signal changes to a low pitch and low cadence this indicates that the DSG® Integration Module for Screw is approaching the cortical wall (Fig.15).


Fig.15

At this time the assembled DSG® Integration Module with the screw may be redirected using anatomical landmarks and discretion. Further drilling without redirection could result in a high pitch and high cadence audio signal that indicates an imminent cortical breach (Fig.16).


Fig.16

Redirect according to the following:

- Change the ratcheting mechanism to neutral position.
- Turn the handle counterclockwise (a high cadence and high pitch audio feedback should be expected, this is due to the blood around the tip of the DSG® Integration Module for Screw).
- Redirect until medium pitch medium cadence signal resumes.
- Engage the ratcheting mechanism in the forward position.
- Continue advancing the device.

Next steps: refer to the partners instructions

4.4 Surgical Technique – MIS or Percutaneous Techniques (if the partner screw system is MIS compatible)

Screw length and diameter selection and connection onto screwdriver: refer to recommendations provided by the partner

The surgeon should locate the intended entry point for the pedicle screws using fluoroscopy, and determine the anticipated angle of entry using his knowledge of MIS procedure for pedicle screw placement including the analysis of pre-op images of the spine of the patient.

Make a short incision on the skin over the targeted pedicle.

Before and during initial insertion of the assembled DSG® Integration Module with the Screw into the pedicle, ensure that the ratcheting mechanism is in the neutral position. This will enhance the ease of penetration into the pedicle.

Remove the Electronic Handle from the DSG® Integration Module for Screw and directly introduce the device through soft tissues. Check the surgical technique of your screw manufacturer for the MIS approach.

Position the assembled DSG® Integration Module with the Screw at the pedicle entry point and along the desired angle of entry using fluoroscopy, and perforate the cortical bone with the tip of the DSG® Integration Module for Screw. Mild malleting on the knob of the DSG® Pin may be required to advance the tip past the entry point.

Then connect the Electronic Handle. As the DSG® Integration Module for Screw is inserted past the entry point into the pedicle a high pitch and high cadence feedback may be heard. This is due to the tip being advanced past the blood collected at the entry point.

Pay close attention to the audio feedback provided as the DSG® Integration Module for Screw is being advanced while applying constant pressure.

Advance the assembled DSG® Integration Module with the Screw into the pedicle along the targeted screw trajectory. Drilling into the pedicle through cancellous bone, a medium pitch and medium cadence audio signal should be heard.

Once a few screw threads have come in contact with bone, engage the ratchet mechanism by turning it clockwise. This will allow to advance the instrument down the pedicle in small increments.

Advance the DSG® integration module for screw assembly through cancellous bone until the screw tip reaches its expected position. Progress of the positioning of the tip of the screw through the body of the vertebral pedicle and into the vertebral body may be routinely checked with fluoroscopy. Muscular contractions may occur if the cortex has been fractured. In this case, the surgeon must stop advancement of the screw, and check and/or alter the initial trajectory.

After each screw insertion, clean the tip of the DSG® Pin with a soft cloth impregnated with saline solution.

In Case Redirection Is Necessary

A change in pitch and cadence of the audio feedback indicates a change in the type of tissue around the tip of the DSG® Integration Module for Screw.

If the audio signal changes to a low pitch and low cadence this indicates that the DSG® Integration Module for Screw is approaching the cortical wall. At this time the assembled DSG® Integration Module with the screw may be redirected using anatomical landmarks and discretion. Further drilling without redirection could result in a high pitch and high cadence audio signal that indicates an imminent cortical breach.

Redirect according to the following:

- Change the ratcheting mechanism to neutral position.
- Turn the handle counterclockwise (a high cadence and high pitch audio feedback should be expected, this is due to the blood around the tip of the DSG® Integration Module for Screw).
- Redirect until medium pitch medium cadence signal resumes.
- Engage ratcheting mechanism in the forward position.
- Continue advancing the device.

Next steps: refer to the partners instructions

5. AFTER USE

The single use instruments must be discarded after the surgical procedure.

Discard the instrument in appropriate containers so that destruction can be performed in a manner that preserves the environment. For detailed information concerning this procedure, please contact your local SpineGuard representative.

All assembled instruments must be disassembled at the point of use before cleaning. Clean all DSG® reusable instruments after each surgery at the point of use with a disposable wipe to remove excessive soil within a maximum of 2 hours from time of soiling.

To disassemble:

1. Remove the Electronic Handle and the DSG® Pin from the DSG® Ratcheting Handle. Dispose of the Electronic Handle and the DSG® Pin (single use only).
2. Separate the screwdriver from the DSG® Ratcheting Handle.

Then follow the “before use”: “cleaning” and “sterilization” procedures mentioned above as soon as possible after use.

If cleaning must be delayed, immerse all DSG® reusable instruments in an enzymatic detergent solution (pH<11) or water at room temperature (15-25°C) to prevent drying of soil and contaminants in and on the device.

6. IMPORTANT MEDICAL INFORMATION

6.1 Indications

The DSG® Integration Module for Screw device is indicated for use during pedicle screw pilot hole drilling to provide feedback to the surgeon via visual and audible alerts that indicate a change in impedance at the tip of the device and may indicate contact of the tip with soft tissues and possible vertebral cortex perforation. DSG® Integration Module for Screw is indicated for use in both open and percutaneous (MIS) surgical approaches to the spine. DSG® Integration Module for Screw is also indicated for use with fluoroscopic guidance in percutaneous (MIS) surgical approaches to the spine. DSG® Integration Module for Screw also is specifically indicated for use in intraoperative electromyographic (“EMG”) surveillance to assist in the location and evaluation of spinal nerves during surgery of the spine, by administration of low voltage electrical energy to tissues and nerves at the operative site, and EMG monitoring of muscle groups associated with those nerves.

6.2 Contraindications

The DSG® Integration Module for Screw is contra-indicated for:

- Pathologies involving the vertebral cortex
- Patients who have received a pacemaker or any other active medical device.

6.3 Warnings

- Do not use the device:
 - In the presence of flammable anesthetics
 - At room temperature above 30°C
 - In a humid environment.
- Use of the DSG® Integration Module for Screw in conditions of extreme osteoporosis is not recommended. The condition of the bone in this situation must be closely evaluated prior to use of the device.
- Only use components provided by SpineGuard® or established partner(s) of SpineGuard. Using components from another manufacturer may cause malfunction of the device and be harmful to the patient.
- The DSG® Pin must be used only with the SpineGuard’s Electronic Handle and the approved partner(s) components.
- No modification of this equipment is allowed
- Intraoperative use of muscle relaxants inhibits muscular contractions which are normally used as a means of detection. Use of muscle relaxants is therefore discouraged when using the DSG® Integration Module for Screw.
- Patient positioning on the operating table is critical: ensure that leg motions (if any) will not cause any injury to the patient or a fall from the operating table.
- Avoid any contact between the DSG® Integration Module for Screw and the patient while using an electrocautery device or defibrillator.
- Use of the device near the thorax may increase the risk of cardiac fibrillation.
- Using the DSG® Integration Module for Screw in association with another electronic device requires special precautions: perturbations on the patient’s ECG monitor may be observed while the DSG® Integration Module for Screw is in use. These abnormalities will stop as soon as the instrument is removed.
- IMPORTANT: Do not allow fluids or foreign bodies in the handle (immersion of the handle, even for a short period of time, is prohibited). In

case of fluids or foreign bodies ingress in the electronics casing, do not use the device.

- Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.
- Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the DSG® Integration Module for Screw. Otherwise degradation of the performance of this equipment could result.
- Operation in close proximity (e.g. 1 m) to a shortwave or microwave therapy ME equipment may produce instability in the device output.

6.4 Precautions

- The choice of the DSG® Pin length and diameter is determined by the pedicle screw system used by the surgeon
- When used as an EMG stimulator, do not let muscular contractions last too long to avoid undue patient fatigue.
- As with any surgical instrument, this device should be carefully inspected before use.
- Always check the integrity of the single use sterile packaging before use. Do not use instrument with opened or damaged packaging.
- If any component of the DSG® Integration Module for Screw is accidentally dropped on the ground during the procedure, do not use it.
- Do not use the DSG® Integration Module for Screw if no rhythmic beep is heard when dipped in a saline solution.
- If the green LED lights are on but do not flash when dipped in a saline solution (circuit is open), do not use the DSG® Integration Module for Screw.
- Always make sure that the "beep" is heard and the green LED lights are on when the tip of the instrument touches the patient wound.
- Do not re-sterilize the DSG® single-use instruments.
- Do not mallet on the Electronic Handle Ref. P2HEXXXX

6.5 Possible adverse effects

The possible adverse effects while using the DSG® Integration Module for Screw are:

- Penetration of the vertebrae resulting in injury or paralysis;
- Muscular contractions induced by nerve root stimulation;
- Loss of neurological functions, appearance of radiculopathies, dural injuries and/or pain. Neurovascular insufficiency, including paralysis or other serious lesions. Cerebrospinal fluid leak;
- Gastrointestinal, urological and/or reproductive tract disorders, including sterility impotence;
- Incapacity to resume the activities of normal everyday life;
- Soft tissue damage;
- Infection;
- Fracture of the vertebral body;
- Death.

7. CLINICAL USE (the clinical study was conducted using the originally cleared PediGuard model (K030526))

A clinical study was performed in which 147 manual pedicle drillings were performed during 28 spinal surgeries using the PediGuard® device in addition to conventional techniques of detecting pedicle screw fracture. Investigators received PediGuard® output when detecting fractures in addition to their conventional surgical methods. A total of 23 vertebral cortex perforations (16%) were confirmed, of which the PediGuard® detected 22. One false positive occurred during the study.

8. STATEMENT (related to European regulation 93/42/EEC)

The DSG® Integration Module for Screw is not intended to be used in direct contact with the Central Nervous System. The contact with any part of the dura mater, may it be the part sheathing the medulla spinalis or the part sheathing the roots of the nerves up to the spinal ganglia, shall be avoided.

9. TECHNICAL SPECIFICATIONS

Storage: Expiry date is mentioned on the single use instruments outer package and must be respected. The DSG® Integration Module for Screw should be stored in a clean dry place.

All the DSG® reusable instruments should be stored in the DSG® Technology Instrument Tray, in a clean dry place.

End of life: If any damage or excessive wear is noted on the DSG® reusable instruments, such as twisting, bending, breakage, pinching at tips, lack of integrity of the etch marks, discoloration or cracks, the DSG® reusable instruments should not be used and should be discarded following the hospital procedures for post-treatment of contaminated devices. In case of complaint, please contact SpineGuard®.

IF YOU HAVE QUESTIONS ABOUT THIS DEVICE, PLEASE CONTACT:

SpineGuard, S.A. (Manufacturer)
10, Cours Louis Lumière
94300 Vincennes France
Phone: +33 (0) 1 45 18 45 19
Fax: +33 (0) 1 45 18 45 20

Symbols found on the products are described below:

See Instructions for Use

Manufacturer

"BF" (body floating) type insulation

Single use only

Expiration date

Sterilized using ethylene oxide

Lot or batch Number

Catalog number

Non-sterile

Information and manufacturer's declaration concerning electromagnetic compatibility (EMC)

Guidance and manufacturer's declaration – electromagnetic emissions		
The PediGuard is intended for use in the electromagnetic environment specified below. The customer or the user of the PediGuard should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment – guidance
RF Emissions CISPR 11	Group 1	PediGuard uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF Emissions CISPR 11	Class A	PediGuard is suitable for use in all establishments other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.

Guidance and manufacturer's declaration – electromagnetic immunity			
The PediGuard is suitable for use in the electromagnetic environment specified below. The customer or the user of the PediGuard should assure that it is used in such an electromagnetic environment.			
Immunity test	Test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 8 kV contact ± 15 kV air	± 8 kV contact ± 15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Radiated RF IEC 61000-4-3	3 V/m from 80 MHz to 2,5 GHz	3 V/m	PediGuard is suitable for use in all establishments other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Radiated RF IEC 61000-4-3 Proximity fields from RF wireless communications equipment	380 - 390 MHz 27 V/m; PM 50%; 18 Hz	27 V/m	Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the PediGuard, including cables specified by the manufacturer.
	430 - 470 MHz 28 V/m; (FM ±5 kHz, 1 kHz sine) PM; 18 Hz	28 V/m	
	704 - 787 MHz 9 V/m; PM 50%; 217 Hz	9 V/m	
	800 - 960 MHz 28 V/m; PM 50%; 18 Hz	28 V/m	
	1700 - 1990 MHz 28 V/m; PM 50%; 217 Hz	28 V/m	
	2400 - 2570 MHz 28 V/m; PM 50%; 217 Hz	28 V/m	
Power frequency magnetic field (50/60Hz) IEC 61000-4-8	5100 - 5800 MHz 9 V/m; PM 50%; 217 Hz	9 V/m	–
	30 A/m	30 A/m	

NOTE The EMISSIONS characteristics of this equipment make it suitable for use in industrial areas and hospitals (CISPR 11 class A). If it is used in a residential environment (for which CISPR 11 class B is normally required) this equipment might not offer adequate protection to radio-frequency communication services. The user might need to take mitigation measures, such as relocating or re-orienting the equipment.

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SpineGuard® S.A.