Instruction For Use of the Straight and Curved PediGuard

LP2-A088-WEB - Version A

Accuracy

Real Time

Safety

Radiation
Free

DSG
Dynamic Surgical Guidance

Educational
Tool

Table of content

<i>EN - User manual</i>	<i>4</i>
<i>FR - Manuel d'utilisation</i>	<i>8</i>
<i>ES - Manual del usuario</i>	<i>12</i>
<i>DE - Benutzerhandbuch</i>	<i>16</i>
<i>IT - Manuale per l'utente</i>	<i>20</i>
<i>PT - Manual do utilizador</i>	<i>24</i>
<i>NL - Handboek voor de gebruiker</i>	<i>28</i>
<i>EL - Οδηγίες χρήσης</i>	<i>32</i>
<i>SK - Návod na použitie</i>	<i>36</i>

USER MANUALEN

This user manual is **not applicable** in the U.S.A. To consult the user manual applicable in the U.S.A. log on to: <https://www.spineguard.com/products-and-applications/spine-products/user-manuals/>

In no case can the PediGuard® replace the surgeon's experience or knowledge of anatomic structures. It is intended to assist the surgeon in the decision-making process during surgery in the operating room.

1. INTRODUCTION

PediGuard is a pedicle preparation device that helps surgeons to drill a pilot hole through the pedicle of the vertebrae during spinal surgery. It serves to alert the surgeons prior to a possible cortex perforation while drilling through the pedicle, by accurately analyzing the electrical conductivity of the surrounding tissues.

2. TECHNICAL DESCRIPTION

The PediGuard is a sterile single use instrument; its handle contains an electronic circuit. The tip of the instrument is the applied part.

The PediGuard is a battery-operated device with a maximum output current of 5.5 milliamps ("mA"). Published literature suggests that threshold level monitoring requires at least 8-10 mA. The PediGuard, which does not perform threshold recordings, is not intended to be used as a substitute for a threshold level monitoring device.

Essential Performances:

1. Mechanics essential performance: the cutting ability during the pedicle drilling.
2. Electronics essential performance: sound and visual signals lead the surgeon during the pedicle drilling. Those informing signals convey the conductivity information which characterizes the tissues in contact with the sensor tip.

DSG Connect only:

An optional feature is available to transmit the conductivity data to an external display to allow signal visualization and recording. The PediGuard devices equipped with this feature are labelled DSG Connect and are equipped with a communication module which transmits data through radio frequency (RF) communication (2.40–2.48 GHz, GFSK modulation and maximum RF output power of 9.9 dBm).

3. FOR USE

3.1 Power-on

To activate the PediGuard, pull the contact tab out of the Handle. Activation is confirmed by a beeping sound and green LEDs.

3.2 Safety features

The PediGuard has several safety features which allow easy checking of proper operation:

Green LEDs	ON	Flashing	OFF	OFF	OFF
Yellow LEDs	OFF	OFF	OFF	OFF	ON
Sound	None	Rated sound	ON	None	None
Case	①	②	③	④	⑤

- ① The instrument is ON. Either the tip of the instrument does not make contact with tissue or is in contact with bone that has very low conductivity.
- ② The tip of the instrument is in contact with tissue.
- ③ The battery is too low: do not use the instrument.
- ④ Either the instrument is not ON, or the battery does not work. Verify that the instrument is ON. If the instrument is ON but the LEDs remain OFF, do not use the instrument.
- ⑤ The instrument tip is not functioning properly: do not use the instrument.

If the LEDs are not visible through the clear wall of the handle, do not use the PediGuard.

3.3 Preliminary checks

Check the PediGuard device for accurate assembly and function by placing the tip of the instrument in saline (0.9% sodium chloride). A high rate high-pitched sound should be heard. The green LEDs should flash at the same frequency. If a metal bowl is used, do not touch the walls of the bowl with the tip of the PediGuard. If no sound is heard, it may be necessary to clean the tip of the instrument using a mild scrubber (e.g. electrocautery scrubber).

3.4 Connection to the DSG Connect App (optional)

The PediGuard devices equipped with the DSG Connect feature can be paired to an external display. The pairing instructions and specific instructions for use are available within the DSG Connect App.

3.5 Surgical technique

Surgical technique – Open Techniques

The surgeon should prepare the intended entry points for the pedicle screws using preferred open techniques to expose the posterior surface of the spine. The correct location for the entry point and angle of entry may be checked by using fluoroscopy. Each hole should be drilled using the following technique:

- Prepare the entry point for drilling by removing the cortical bone.
- Direct the tip of the PediGuard to the soft tissue (wound) to initiate a rated high-pitched "beep".
- Position the PediGuard at the entry point and advance its tip through cancellous bone in the expected position of the screw. Use constant hand pressure and gentle alternate rotating motions (*never use mallet blows*).

Signal analysis

A change in pitch and cadence of the audio feedback and LEDs flashing frequency indicates a change in conductivity around the tip of the PediGuard. Each type of biologic tissue presents a specific electrical conductivity. Blood presents a high electrical conductivity at the opposite of cortical bone that presents a low conductivity, while cancellous bone presents a medium electrical conductivity. Drilling through cancellous bone, a medium pitched beep should be heard. If the audio signal changes to a low rated pitched beep, this may indicate that the device is approaching or entering the cortical wall. If the audio signal changes to a high rated pitched beep, it may indicates a bone breach.

Muscular contractions may occur if the cortex has been fractured. In this case, the surgeon must stop advancement of the tip, and check the trajectory. If redirection is necessary, instrument should be pulled back before to adjust trajectory.

Before drilling a new hole, clean the tip of the instrument with a soft cloth impregnated with saline solution.

3.6 After surgery

The instrument must be discarded after the surgical procedure.

Discard the instrument in appropriate containers so that destruction can be performed in a manner that preserves the environment. For detailed information concerning this procedure, please contact your local SpineGuard representative.

4. IMPORTANT MEDICAL INFORMATION

4.1 Intended use (intended purpose)

The PediGuard is intended to drill pilot holes for pedicle screws. It provides feedback to the surgeon via visual and audible signal that indicate a change in conductivity at the tip of the probe and may indicate contact of the tip with soft tissues and possible vertebral cortex perforation.

4.2 Indications

The PediGuard is indicated for all surgical procedures requiring the drilling of pilot holes for pedicle screws.

The PediGuard is indicated for use in open surgical approaches to the spine.

4.3 Intended patient population

The PediGuard is intended for skeletally mature patient, either adult or pediatric, undergoing spine surgery requiring the drilling of vertebrae for the insertion of pedicle screws. The PediGuard Gearshift cannot be used on premature or neonatal infants, and patients must have a minimum weight of 25 kg.

4.4 Intended environment (users)

The PediGuard is intended to be used by surgeons in hospitals.

4.5 Contraindications

- Pathologies involving the vertebral cortex.
- Patients who have received a pacemaker or any other active medical device.
- The PediGuard is not intended to be used in direct contact with the central nervous system.
- The contact with any part of the dura mater, may it be the part sheathing the medulla spinalis or the part sheathing the roots of the nerves up to the spinal ganglia, shall be avoided.

4.6 Warnings

- Do not use the device:
 - In the presence of flammable anesthetics
 - At room temperature above 30°C
 - In a humid environment.
- The use of the PediGuard in conditions of extreme osteoporosis is not recommended. The condition of the bone in this situation must be closely evaluated prior to use of the device.
- Only use components provided by SpineGuard. Using components from another manufacturer may cause malfunction of the device and be harmful to the patient.
- No modification of this equipment is allowed.
- Patient positioning on the operating table is critical: ensure that leg motions (if any) will not cause any injury to the patient or a fall from the operating table.
- Avoid any contact between the PediGuard and the patient while using an electrocautery device or defibrillator.
- Use of the device near the thorax may increase the risk of cardiac fibrillation.
- Using the PediGuard in association with another electronic device requires special precautions: perturbations on the patient's ECG monitor may be observed while the PediGuard is in use. These abnormalities will stop as soon as the instrument is removed.
- IMPORTANT: Do not allow fluids or foreign bodies in the handle (immersion of the handle, even for a short period of time, is prohibited). In case of fluids or foreign bodies ingress in the electronics casing, do not use the device.
- Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.
- Portable RF communications equipment (including peripherals such as antenna cables and externals antennas) should be used no closer than 30 cm (12 inches) to any part of the PediGuard System. Otherwise degradation of the performance of this equipment could result.
- Operation in close proximity (e.g. 1m) to a shortwave or microwave therapy ME equipment may produce instability in the device output.

4.7 Precautions

- The choice of the PediGuard diameter is determined by the diameter of the screws the surgeon wants to place and by the pedicle size.
- As with any surgical instrument, this device should be carefully inspected before use.
- Before opening, always perform a visual check of the integrity of the sterile pouch. It must not show any damages or perforations and must be sealed. If the pouch is opened or damaged, the product must not be used.
- If the PediGuard is accidentally dropped on the ground, do not use it.
- Do not use the PediGuard if no rhythmic beep is heard when dipped in a saline solution.
- If the green LED lights are ON but do not flash when dipped in a saline solution (circuit is open), do not use the PediGuard.
- Always make sure that the "beep" is heard and the green LED lights are on when the tip of the instrument touches the patient wound.
- Do not re-sterilize the instrument.

4.8 Possible adverse effects

The possible adverse effects while using the PediGuard system are:

Adverse effects identified in clinical studies which represent an occurrence rate of less than 0.5%:

- Neurological complications: loss of neurological functions, appearance of radiculopathies, dural injuries and/or pain. Neurovascular insufficiency, including paralysis or other serious lesions,

Adverse effects identified in clinical studies which represent an occurrence rate of 0% in the scientific literature and post-market surveillance:

- Biomechanical complications: muscular contractions induced by nerve root stimulation; fracture of the vertebral body,
- Cardiological complications, particularly in patients with pacemakers,
- Cerebrospinal fluid leak,
- Vascular complications.
- Other (gastrointestinal, urological and/or reproductive tract disorders, infection, soft tissue damage).

5. TECHNICAL SPECIFICATIONS

Storage: Expiry date is mentioned on the outer package and must be respected. The PediGuard should be stored in a clean dry place and in ambient conditions.

Service life: Once ON, the PediGuard cannot be turned OFF. The PediGuard battery provides sufficient autonomy.

6. INFORMATION RELATED TO SINGLE USE DEVICES

PediGuard is a single use device. The retreatment and the re-use of the PediGuard are forbidden.

Several potential dangers related to the re-use of the PediGuard can increase the risk to the patient. Some technical characteristics of the PediGuard are not compatible with the re-use of the instrument. Possible issues with re-sterilization and re-use are: non-functioning device, failure of some or all of the electronic functions, unknown changes to the electronic function of the device, non-sterile device, transfer of disease, cleaning agents or disinfectants.

Additionally, the re-use of a single use device will result in the loss of traceability of the medical device and the loss of the technical documentation, as the instruction for use.

Important information: any serious incident relating to the device must be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is established. For such a report, or for any other problem or malfunction, please contact: complaints@spineguard.com

ONLY USE INSTRUMENTS WITH INTACT PACKAGING. DO NOT CLEAN, RE-STERILIZE, OR RE-USE THIS DEVICE. IF YOU HAVE QUESTIONS ABOUT THIS DEVICE, PLEASE CONTACT:

SpineGuard, S.A. (MANUFACTURER)
10, Cours Louis Lumière
94300 Vincennes - France
Phone: +33 (0) 1 45 18 45 19
Fax: +33 (0) 1 45 18 45 20

Symbols found on the product are described below:

See instructions for use

Manufacturer

Single use only

Refer to instruction manual

“BF” (body floating) type insulation

Expiration date

Double sterile barrier system

Medical device

Lot or batch Number

Catalog number

Information and manufacturer's declaration concerning electromagnetic compatibility (EMC)

Guidance and manufacturer's declaration – electromagnetic emissions		
The PediGuard is intended for use in the electromagnetic environment specified below. The customer or the user of the PediGuard should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment – guidance
RF Emissions CISPR 11	Group 1	PediGuard uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF Emissions CISPR 11	Class A	PediGuard is suitable for use in all establishments other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.

Guidance and manufacturer's declaration – electromagnetic immunity			
The PediGuard is suitable for use in the electromagnetic environment specified below. The customer or the user of the PediGuard should assure that it is used in such an electromagnetic environment.			
Immunity test	Test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 8 kV contact ± 15 kV air	± 8 kV contact ± 15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Radiated RF IEC 61000-4-3	3 V/m from 80 MHz to 2,5 GHz	3 V/m	PediGuard is suitable for use in all establishments other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.

Guidance and manufacturer's declaration – electromagnetic immunity			
The PediGuard is suitable for use in the electromagnetic environment specified below. The customer or the user of the PediGuard should assure that it is used in such an electromagnetic environment.			
Immunity test	Test level	Compliance level	Electromagnetic environment – guidance
Radiated RF IEC 61000-4-3 Proximity fields from RF wireless communications equipment	380 - 390 MHz 27 V/m; PM 50%; 18 Hz	27 V/m	Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the PediGuard, including cables specified by the manufacturer.
	430 - 470 MHz 28 V/m; (FM ±5 kHz, 1 kHz sine) PM; 18 Hz	28 V/m	
	704 - 787 MHz 9 V/m; PM 50%; 217 Hz	9 V/m	
	800 - 960 MHz 28 V/m; PM 50%; 18 Hz	28 V/m	
	1700 - 1990 MHz 28 V/m; PM 50%; 217 Hz	28 V/m	
	2400 - 2570 MHz 28 V/m; PM 50%; 217 Hz	28 V/m	
Power frequency magnetic field (50/60Hz) IEC 61000-4-8	5100 - 5800 MHz 9 V/m; PM 50%; 217 Hz	9 V/m	
	30 A/m	30 V/m	

The EMISSIONS characteristics of this equipment make it suitable for use in industrial areas and hospitals (CISPR 11 class A). If it is used in a residential environment (for which CISPR 11 class B is normally required) this equipment might not offer adequate protection to radio-frequency communication services. The user might need to take mitigation measures, such as relocating or re-orienting the equipment.

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SpineGuard® S.A.