Instruction For Use of the Electronic Handle

LP2-A110-Web - Version A

Accuracy

Real Time

Safety

Educational
Tool

Radiation
Free

DSG
Dynamic Surgical Guidance

Table of content

<i>EN - User manual</i>	4
<i>FR - Manuel d'utilisation</i>	9
<i>ES - Manual del usuario</i>	15
<i>DE - Benutzerhandbuch</i>	21
<i>IT - Manuale per l'utente</i>	27
<i>PT - Manual do utilizador</i>	32
<i>NL - Handboek voor de gebruiker</i>	37
<i>EL - Οδηγίες χρήσης</i>	43
<i>SK - Návod na použitie</i>	49

USER MANUAL

EN

This user manual is **not applicable** in the U.S.A. To consult the user manual applicable in the U.S.A. log on to: <https://www.spineguard.com/products-and-applications/spine-products/user-manuals/>

The Electronic Handle is a PediGuard® medical device integrated into the following configurable medical devices: the Cannulated PediGuard and the PediGuard Threaded. In no case can the PediGuard replace the surgeon's experience or knowledge of anatomic structures. It is intended to assist the surgeon in the decision-making process during surgery in the operating room.

The Electronic Handle is a medical device available in two versions:

- The DSG Handle (reference P2HE1000)
- The DSG Connect Handle (reference P2HE2000)

Unless otherwise indicated, this user manual applies to both references which will be identified as: the Electronic Handle.

For use with the Cannulated PediGuard, part A and part B of this user manual apply.

For use with the Threaded PediGuard, part A and part C of this user manual apply.

Part A : The Electronic Handle

1. GENERAL DESCRIPTION

The Electronic Handle is an electrical device designed to be used in surgery in combination with SpineGuard probes in drilling pilot holes into pedicles and vertebral bodies. It serves to alert the surgeon to a possible vertebral cortex fracture while drilling a pedicle screw pilot hole.

2. TECHNICAL DESCRIPTION

The Electronic Handle contains an electronic circuit and comes in the form of a ready-to-use sterile single-use instrument.

The Electronic Handle is a battery-operated device with a maximum output current of 5.5 milliamps ("mA"). Published literature suggests that threshold level monitoring requires at least 8-10 mA. The Electronic Handle, which does not perform threshold recordings, is not intended to be used as a substitute for a threshold level monitoring device.

Essential Performances:

Electronics essential performances: Sound and visual signals lead the surgeon during the pedicle drilling. Those informing signals convey the conductivity information which characterizes the tissues in contact with the sensor tip.

About the DSG Connect Handle:

An optional feature is available to transmit conductivity data to an external display to allow signal visualization and recording. The devices equipped with this feature are labelled DSG Connect and are equipped with a communication module which transmits data through radio frequency (RF) communication (2.40–2.48 GHz, GFSK modulation and maximum RF output power of 9.9 dBm).

3. INSTRUCTION FOR USE OF THE ELECTRONIC HANDLE

3.1 Assembly instruction

Once you have removed the Electronic Handle from its sterile packaging, it must be assembled with the other configurable devices making up the PediGuard Cannulated (refer to Part B) or the PediGuard Threaded (refer to part C).

3.2 Power-on

To activate the Electronic Handle, pull the activation tab out of the Handle. Activation is confirmed by a beeping sound and green LEDs.

3.3 Safety features

Once assembled with the other devices making up the Cannulated PediGuard or the Threaded PediGuard, the Electronic Handle is fitted with safety devices to check that it is working properly.

The PediGuard has several safety features which allow easy checking of proper operation:

Green LEDs	ON	Flashing	OFF	OFF	OFF
Yellow LEDs	OFF	OFF	OFF	OFF	ON
Sound	None	Rated sound	Continuous crackling	None	None
Case	①	②	③	④	⑤

① The instrument is ON. Either the tip of the instrument does not make contact with tissue or is in contact with bone that has very low conductivity.

② The tip of the instrument is in contact with tissue.

③ The battery is too low: do not use the instrument.

④ Either the instrument is OFF, or the battery does not work. Verify that the activation tab is removed. If the instrument is ON but the LEDs remain OFF, do not use the instrument.

⑤ The instrument tip is not functioning properly: do not use the instrument.

If the LEDs are not visible through the clear wall of the handle, do not use the Electronic Handle.

3.4 Preliminary checks

Once the electronic handle has been assembled with the desired configurable devices, place the tip of the instrument in saline (0.9% sodium chloride). A high rate high-pitched sound should be heard. The green LEDs should flash at the same frequency. If a metal bowl is used, do not touch the bowl.

If no sound is heard, it may be necessary to clean the tip of the instrument using a mild scrubber (e.g. electrocautery scrubber) or gently tap the handle of the PediGuard in the palm of your hand.

3.5 Connection to the DSG Connect App (optional) only applies to DSG Connect Handle

The DSG Connect Handle equipped with the DSG Connect feature can be paired to an external display. The pairing instructions and specific instructions for use are available within the DSG Connect App.

3.6 After surgery

The electronic handle must be discarded after the surgical procedure. Discard the instrument in appropriate containers so that destruction can be performed in a manner that preserves the environment. For detailed information concerning this procedure, please contact your local SpineGuard representative.

Important note: in Chapters 4, 5 and 6 the «Device» refers to the electronic handle assembled with other dedicated medical devices in Cannulated PediGuard or Threaded PediGuard configuration.

4. IMPORTANT MEDICAL INFORMATION

4.1 Intended Use

The Device is intended for use during pedicle screw pilot hole drilling to provide feedback to the surgeon via visual and audible alerts that indicate a change in impedance at the tip of the probe and may indicate contact of the tip with soft tissues and possible vertebral cortex perforation.

The Threaded PediGuard is also intended during pilot-holes drilling to receive screws in the vertebral body (anterior surgical approach).

4.2 Indications

The Device is indicated for all surgical procedures requiring the drilling of pilot holes for pedicle screws.

The Threaded PediGuard is also indicated for drilling pilot-holes to receive screws in the vertebral body.

The Device is indicated for use in both open and percutaneous (MIS) surgical approaches to the spine. The Device is also indicated for use with fluoroscopic guidance in percutaneous (MIS) surgical approaches to the spine.

4.3 Targeted Population

The Device is intended for adult and pediatric patients undergoing spinal surgery requiring the drilling of pilot holes in the vertebrae. The Device cannot be used on premature or newborn babies, and patients must weigh at least 25 kg.

4.4 Intended Environment (users)

The Device is intended to be used by surgeons in hospitals.

4.5 Contraindications

- Pathologies involving the vertebral cortex
- Patients who have received a pacemaker or any other active medical device.
- The Device must not be used in direct contact with the central nervous system.
- The contact with any part of the dura mater, may it be the part sheathing the medulla spinalis or the part sheathing the roots of the nerves up to the spinal ganglia, shall be avoided.

4.6 Warnings

- Do not use the device:
 - In the presence of flammable anesthetics,
 - At room temperature above 30°C,
 - In a humid environment.
- Use of the Device in conditions of extreme osteoporosis is not recommended. The condition of the bone in this situation must be closely evaluated prior to use of the device.
- Only use components provided by SpineGuard. Using components from another manufacturer may cause malfunction of the device and may be harmful to the patient.
- No modification of this equipment is allowed.
- Patient positioning on the operating table is critical: ensure that leg motions (if any) will not cause any injury to the patient or a fall from the operating table.
- Avoid any contact between the Device and the patient while using an electrocautery device or defibrillator.
- Use of the device near the thorax may increase the risk of cardiac fibrillation.
- Using the Device in association with another electronic device requires special precautions: perturbations on the patient's ECG monitor may be observed while the Device is in use. These abnormalities will stop as soon as the instrument is removed.
- IMPORTANT: Do not allow fluids or foreign bodies in the electronic handle (immersion of the electronic handle, even for a short period of time, is prohibited). In case of fluids or foreign bodies ingress in the electronics casing, do not use the device.
- Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.
- Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the Device. Otherwise degradation of the performance of this equipment could result.
- Operation in close proximity (e.g. 1m) to a shortwave or microwave therapy ME equipment may produce instability in the device output.

4.7 Precautions

- The choice of the Device diameter and configuration is determined by the diameter of the screws the surgeon wants, the anatomy to place (pedicle or vertebral body) and by the size of the vertebrae.
- As with any surgical instrument, this device should be carefully inspected before use.
- Before opening, always perform a visual check of the integrity of the sterile pouch. It must not show any damages or perforations and must be sealed. If the pouch is opened or damaged, the product must not be used.
- If the PediGuard is accidentally dropped on the ground, do not use it.
- Do not use the PediGuard if no rhythmic beep is heard when dipped in a saline solution.
- If the green LED lights are ON but do not flash when dipped in a saline solution (circuit is open), do not use the PediGuard.
- Always make sure that the "beep" is heard and the green LED lights are on when the tip of the instrument touches the patient wound.
- Do not re-sterilize the Electronic Handle.
- For the anterior approach, the lungs are selectively ventilated allowing the lung on the convex side to collapse creating the workspace, however there is a risk of perforation of lungs if not properly collapsed.
- Light impacts are acceptable on the Electronic Handle, but it is recommended not to tap excessively in order to avoid damaging the electronic system.
- Before drilling a new hole, clean the tip of the instrument with a soft cloth impregnated with saline solution.

Note: whatever the surgical technique used, or the configuration of the medical device used, it is possible that muscular contractions may occur if the bone cortex is breached. In this case, stop the progression of the instrument and check the trajectory.

4.8 Possible Adverse Effects

The possible adverse effects while using the PediGuard system are:

Adverse effects identified in clinical studies which represent an occurrence rate of less than 0.5%:

- Neurological complications: loss of neurological functions, appearance of radiculopathies, dural injuries and/or pain. Neurovascular insufficiency, including paralysis or other serious lesions,
- Other (gastrointestinal, urological and/or reproductive tract disorders, infection, soft tissue damage).

Adverse effects identified in clinical studies which represent an occurrence rate of 0% in the scientific literature and post-market surveillance:

- Biomechanical complications: muscular contractions induced by nerve root stimulation; fracture of the vertebral body,
- Cardiological complications, particularly in patients with pacemakers,
- Cerebrospinal fluid leak,
- Vascular complications.

5. TECHNICAL SPECIFICATIONS

Storage: Expiry date is mentioned on the outer package and must be respected. The Electronic Handle should be stored in a clean dry place under ambient temperature conditions.

Service life: Once ON, the Electronic Handle cannot be turned OFF. The Electronic Handle battery provides sufficient autonomy for a standard spine surgery.

6. INFORMATION RELATED TO SINGLE USE DEVICES

The Electronic Handle is a single use device. The retreatment and the re-use of the Electronic Handle are forbidden.

Several potential dangers related to the re-use of the Electronic Handle can increase the risk to the patient. Some technical characteristics of the Electronic Handle are not compatible with the re-use of the instrument. Possible issues with re-sterilization and re-use are:

- Non-functioning device,
- Failure of some or all of the electronic functions,
- Unknown changes to the electronic function of the device,
- Non-sterile device,
- Transfer of disease, cleaning agents or disinfectants.

Additionally, the re-use of a single use device will result in the loss of traceability of the medical device and the loss of the technical documentation, as the instruction for use.

Symbols found on the product are described below:

**Information and manufacturer's declaration concerning electromagnetic compatibility (EMC)**

Guidance and manufacturer's declaration – electromagnetic emissions		
The Device is intended for use in the electromagnetic environment specified below. The customer or the user of the Device should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment – guidance
RF Emissions CISPR 11	Group 1	The Device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF Emissions CISPR 11	Class A	The Device is suitable for use in all establishments other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.

Guidance and manufacturer's declaration – electromagnetic immunity			
The Device is suitable for use in the electromagnetic environment specified below. The customer or the user of the Device should assure that it is used in such an electromagnetic environment.			
Immunity test	Test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 8 kV contact ± 15 kV air	± 8 kV contact ± 15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Radiated RF IEC 61000-4-3	3 V/m from 80 MHz to 2,5 GHz	3 V/m	The Device is suitable for use in all establishments other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Radiated RF IEC 61000-4-3 Proximity fields from RF wireless communications equipment	380 - 390 MHz 27 V/m; PM 50%; 18 Hz	27 V/m	Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the PediGuard, including cables specified by the manufacturer.
	430 - 470 MHz 28 V/m; (FM ±5 kHz, 1 kHz sine) PM; 18 Hz	28 V/m	
	704 - 787 MHz 9 V/m; PM 50%; 217 Hz	9 V/m	
	800 - 960 MHz 28 V/m; PM 50%; 18 Hz	28 V/m	
	1700 - 1990 MHz 28 V/m; PM 50%; 217 Hz	28 V/m	
	2400 - 2570 MHz 28 V/m; PM 50%; 217 Hz	28 V/m	
	5100 - 5800 MHz 9 V/m; PM 50%; 217 Hz	9 V/m	

Guidance and manufacturer's declaration – electromagnetic immunity			
The Device is suitable for use in the electromagnetic environment specified below. The customer or the user of the Device should assure that it is used in such an electromagnetic environment.			
Immunity test	Test level	Compliance level	Electromagnetic environment – guidance
Power frequency magnetic field (50/60Hz) IEC 61000-4-8	30 A/m	30 V/m	-

The EMISSIONS characteristics of this equipment make it suitable for use in industrial areas and hospitals (CISPR 11 class A). If it is used in a residential environment (for which CISPR 11 class B is normally required) this equipment might not offer adequate protection to radio-frequency communication services. The user might need to take mitigation measures, such as relocating or re-orienting the equipment.

Part B: The Electronic Handle used in the configurable Cannulated PediGuard device

In this configuration, the Electronic Handle must be assembled with the P2NDXXXX Needle, a single-use device manufactured by SpineGuard.

It cannot be assembled with other medical devices, manufactured by other manufacturers than SpineGuard.

Assembly instruction

- Before connecting the P2NDXXXX Needle to the handle, check that the active stylet (orange cap) is fully screwed into the cannula. The distal tip of the active stylet should be visible at the end of the cannula.
- Connect the Needle into the Electronic Handle. You should feel it click into place. Check that the needle is fully engaged in the Electronic Handle.
- Activate the device by pulling the contact tab out of the Electronic Handle. Activation is confirmed by a beeping sound and green LEDs lighting up.
- The Cannulated PediGuard is ready for use.

Surgical technique – Open Techniques

The surgeon should prepare the intended entry points for the pedicle screws using preferred open techniques to expose the posterior surface of the spine. The correct location for the entry point and angle of entry may be checked by using fluoroscopy.

Each hole should be drilled using the following technique:

- Prepare the entry point for drilling by removing the cortical bone.
- Direct the tip of the Cannulated PediGuard to the soft tissue (wound) to initiate a rated high-pitched "beep".
- Position the Cannulated PediGuard at the entry point and advance its tip through cancellous bone in the expected position of the screw. Use constant hand pressure and gentle alternate rotating motions. A medium pitch, medium cadence sound must be heard.

Surgical technique – MIS or percutaneous Techniques

The surgeon should locate the intended entry point for the pedicle screws using fluoroscopy, and determine the anticipated angle of entry using his knowledge of MIS procedure for pedicle screw placement including the analysis of pre-op images of the spine of the patient.

Each hole should be drilled using the following technique:

- In the case of a percutaneous procedure, directly introduce the Cannulated PediGuard through soft tissues.
- In the case of a mini-open procedure, use the specific retractor provided with the pedicle screw instrumentation set until the entry point can be directly visualized.
- Position the Cannulated PediGuard at the pedicle entry point and along the desired angle of entry using fluoroscopy, and perforate the cortical bone with the tip of the Cannulated PediGuard.
- Advance the Cannulated PediGuard tip through cancellous bone in the expected position of the screw. Use constant hand pressure and gentle alternate rotating motions. A medium pitch, medium cadence sound must be heard.
- Progress of the positioning of the tip of the Cannulated PediGuard through the body of the pedicle and into the vertebral body may be routinely checked with fluoroscopy.
- A change in pitch and cadence of the audio feedback indicates a change in the type of tissue around the tip of the Cannulated PediGuard. If the audio signal changes to a low pitch and low cadence this indicates that the Cannulated PediGuard is approaching the cortical wall. At this time the Cannulated PediGuard may be redirected using anatomical landmarks and discretion.

Part C: The Electronic Handle used in the configurable Threaded PediGuard device

In this configuration the Electronic Handle must be assembled with the following devices manufactured exclusively by SpineGuard:

Single-use instruments provided sterile: DSG® Pin Ref. D1PUXXXX

Reusable instruments that must be cleaned and sterilized before use:

- DSG® Ratcheting Handle – Ref. D1RHXXXX
- DSG® Threaded Shafts – Ref. D1TAXXXX
- DSG® Sleeve – Ref. D1STXXXX
- DSG® Technology Instrument Tray: Base Ref. D1TT1XXX and Lid Ref. D1TT2XXX.

It cannot be assembled with other medical devices manufactured by manufacturers other than SpineGuard.

Assembly instructions

- Attach the selected Threaded Shaft to the Ratcheting Handle. Check that the Threaded Shaft is fully engaged with the Ratcheting Handle - the alignment mark on the Threaded Shaft should be flush with the distal end of the Ratcheting Handle.
- Insert the Pin through the Ratcheting Handle - Threaded Shaft assembly. Check that the knob of the Pin is fully screwed into the Ratcheting Handle. The distal tip of the Pin should be visible (3mm +/- 1mm) at the distal end of the Threaded Shaft.
- Connect the Electronic Handle to the DSG Pin - Threaded Shaft assembly. Ensure that the Electronic Handle is securely connected to the Ratcheting Handle.
- Activate the Threaded PediGuard device by pulling the contact tab out of the Electronic Handle. Activation is confirmed by a beeping sound and green LEDs lighting up.
- In the case of a percutaneous procedure, the DSG Sleeve may be used with the Threaded PediGuard device by inserting the Threaded Shaft (with pin) through the sleeve.
- The Threaded PediGuard is now ready for use.

For further information on assembly instructions, please refer to the assembly instructions provided in the DSG Pin box or on the SpineGuard website, section «User manuals: Threaded PediGuard».

Surgical techniques– Open and minimally invasive techniques

The surgeon should prepare the intended entry points for the pedicle screws using preferred open techniques to expose the posterior surface of the spine. The correct location for the entry point and angle of entry may be checked by using fluoroscopy. For each drilling, scrupulously follow the procedure below:

- In the case of a minimally invasive approach, use the retractor provided with the screw instrumentation set until the entry point can be directly visualized.
- Prepare the entry point for pedicle drilling by removing the cortical bone using conventional techniques.
- Before and during the initial insertion of the PediGuard Threaded into the bone, ensure that the ratcheting mechanism is in the neutral position.
- Position the PediGuard Threaded at the entry point. A high-pitched and high cadence sound should be heard.
- Advance the PediGuard Threaded into the pedicle along the targeted screw trajectory. Drilling into the pedicle through cancellous bone, a medium pitch and medium cadence audio signal should be heard
- Once the first threads of the PediGuard Threaded are engaged in the bone, engage the ratcheting handle mechanism by turning it clockwise. This will allow to advance the instrument down the pedicle in small increments through a series of alternating rotational movements of the handle.
- The graduated marks on the threaded shaft indicate the progression of its distal tip into the pedicle and can be used to estimate the length of screws to be inserted. Note that the etch marks on the Threaded Shaft include ~3mm for the tip of the pin at the distal end of the threaded shaft.
- A change in pitch and cadence of the audio feedback indicates a change in the type of tissue around the tip of the PediGuard Threaded device. If the audio signal changes to a low pitch and low cadence this indicates that the PediGuard Threaded device is approaching the cortical wall. At this time, the PediGuard Threaded may be redirected using anatomical landmarks. Redirect according to the following:
 - Change the ratcheting mechanism to neutral position.
 - Turn the PediGuard Threaded counterclockwise. A high cadence and high pitch audio feedback should be expected, this is due to the blood around the tip.
 - Redirect until medium pitch medium cadence signal resumes.
 - Engage ratcheting mechanism in the forward position and continue advancing the PediGuard Threaded.

Surgical Technique – Percutaneous Techniques

The surgeon should locate the intended entry point for the pedicle screws using fluoroscopy, and determine the anticipated angle of entry using his knowledge of MIS procedure for pedicle screw placement including the analysis of pre-op images of the spine of the patient.

For each drilling, scrupulously follow the procedure below:

- Before and during initial insertion of the assembly into the pedicle, ensure that the ratcheting mechanism is in the neutral position.
- Remove the Electronic Handle from the Threaded PediGuard device and directly introduce the Threaded PediGuard device through soft tissues. The DSG Sleeve may be used to protect the soft tissues: load the DSG Sleeve onto the DSG Threaded Shaft before introducing the Threaded PediGuard device through soft tissues. Note: Do not use the DSG Sleeve with DSG Threaded Shafts with diameters larger than 7.0mm.
- Position the Threaded PediGuard device at the pedicle entry point and along the desired angle of entry using fluoroscopy, and perforate the cortical bone with the tip of the Threaded PediGuard device. Mild malleting on the knob (of the DSG Pin) may be required to advance the tip past the entry point.
- When the DSG Sleeve is used with the Threaded PediGuard device:
 - Slide down the sleeve to put it in contact with the bone.
 - The graduated marks on the threaded shaft, visible above the Sleeve, indicate the progression of the distal tip of the instrument into the pedicle. The double lines on the threaded shaft indicate that the first thread has been engaged with the bone enabling the use of the ratcheting mechanism
 - Then connect the Electronic Handle. As the Threaded PediGuard device is inserted past the entry point into the pedicle a high pitch and high cadence feedback may be heard
- Advance the Threaded PediGuard device into the pedicle along the targeted screw trajectory. Drilling into the pedicle through cancellous bone, a medium pitch and medium cadence audio signal should be heard.
- After the tip of the Threaded PediGuard device has been inserted into the pedicle and the first thread has engaged with the bone, enable the ratcheting mechanism by turning it clockwise. This will allow to advance the instrument down the pedicle in small increments.
- Advance its tip through cancellous bone in the expected position of the screw. The marks on the DSG® Threaded Shaft can be used to estimate the length of the screw to be inserted. Note that the etch marks on the Threaded Shaft include ~3mm for the tip of the pin at the distal end of the Threaded PediGuard device. When the DSG Sleeve is being used, the etch marks on the Threaded Shaft, visible above the Sleeve, should be read. Progress of the positioning of the tip of the Threaded PediGuard device through the body of the vertebral pedicle and into the vertebral body may be routinely checked with fluoroscopy.
- A change in pitch and cadence of the audio feedback indicates a change in the type of tissue around the tip of the Threaded PediGuard device. If the audio signal changes to a low pitch and low cadence this indicates that the Threaded PediGuard device is approaching the cortical wall. At this time the Threaded PediGuard device may be redirected using anatomical landmarks and discretion. Redirect according to the following:
 - Change the ratcheting mechanism to neutral position.
 - Turn the Threaded PediGuard device counterclockwise (a high cadence and high pitch audio feedback should be expected – this is due to the blood around the tip).
 - Redirect until medium pitch medium cadence signal resumes.
 - Engage ratcheting mechanism.
 - Continue advancing the Threaded PediGuard device.
- Once desired trajectory and depth is reached detach the electronic handle. The DSG Pin is then removed to insert a guidewire through the threaded shaft – ratcheting handle assembly into the pedicle. Be extremely careful with the position of the guidewire. Unintentional advancement of the wire can potentially be very dangerous. Once the guidewire is inserted, reverse the direction of the ratcheting mechanism and unscrew the Threaded Shaft.

Surgical Technique – Vertebral body pilot hole anterior approach

- Prior to drilling the vertebrae with the Threaded PediGuard device, prepare the entry point and perforate the vertebral cortex using standard techniques.
- Before and during initial insertion of the assembly into the vertebrae, ensure that the ratcheting mechanism is in the neutral position. This will enhance the ease of penetration into the vertebrae.
- Direct the tip of the Threaded PediGuard device to the soft tissue (wound) to initiate a rated high-pitched “beep”. The DSG Sleeve may be used to protect soft tissue.
- After the tip of the Threaded PediGuard device has been inserted into the vertebrae and the first thread has engaged with the bone, enable the ratcheting mechanism by turning it clockwise. This will allow to advance the instrument down the vertebrae in small increments.
- Advance the tip through cancellous bone in the expected position of the screw. The marks on the DSG® Threaded Shaft can be used to estimate the length of the screw to be inserted. Note that the etch marks on the Threaded Shaft include ~3mm for the tip of the pin at the distal end of the Threaded PediGuard device.
- A change in pitch and cadence of the audio feedback indicates a change in the type of tissue around the tip of the Threaded PediGuard device. If the audio signal changes to a low pitch and low cadence this indicates that the Threaded PediGuard device is approaching the cortical wall. If deemed necessary, the Threaded PediGuard device may be redirected using anatomical landmarks and discretion. Redirect according to the following:
 - Change the ratcheting mechanism to neutral position.
 - Turn the Threaded PediGuard device counterclockwise (a high cadence and high pitch audio feedback should be expected – this is due to the blood around the tip).
 - Redirect until medium pitch medium cadence signal resumes.
 - Engage ratcheting mechanism.
 - Continue advancing the Threaded PediGuard device.

Important information: Any serious incident relating to the device must be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is established. For such a report, or for any other problem or malfunction, please contact: complaints@spineguard.com

ONLY USE INSTRUMENTS WITH INTACT PACKAGING. DO NOT CLEAN, RE-STERILIZE, OR RE-USE THIS DEVICE. IF YOU HAVE QUESTIONS ABOUT THIS DEVICE, PLEASE CONTACT:

SpineGuard, S.A. (Manufacturer)
10, Cours Louis Lumière
94300 Vincennes - France
Phone: +33 (0) 1 45 18 45 19
Fax: +33 (0) 1 45 18 45 201.

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SpineGuard, S.A. (Výrobca)
SpineGuard® S.A.